## Table of Contents

<b>ENGLISH</b> .....	<b>2</b>
.....	<b>8</b>

This document is an electronic copy of the basic 'Instructions for Use' of Simpleware Medical Software. Detailed instructions can be found in the 'Reference Guide'.

## 1 General Description of Simpleware Medical

Simpleware™ Medical is standalone software package for converting 3D image data into models that can be exported in different file formats for applications as per the intended use.

## 2 Intended Use

Simpleware Medical is intended for use as a software interface and image segmentation system for the transfer of medical imaging information to an output file. The software is intended to inform clinical management by improving confidence in clinical decision-making and pre-surgical planning through use of patient-specific 3D models (virtual 3D models and 3D printed anatomical models) created from medical imaging. Outputs from the software must be reviewed for clinical suitability by expert clinicians working on specific patient cases.

### 2.1 Clinical Benefits

Simpleware Medical provides potential indirect benefits to intended users:

- Improved understanding of individual patient anatomies through 3D model visualization and segmentation as an improvement on 2D visualization
- Improved planning of surgeries through virtual 3D models and 3D printed anatomical models
- Improved ability to select safest/optimal treatment (implant choice/device sizing and customization to specific patients) through accurate measurements of anatomy
- Improved future clinical procedures through analysis of past procedures by comparing pre-surgical 3D

## 2.3 Patient Population

Simpleware Medical has no intended patient population. Outputs (virtual 3D models and 3D printed anatomical models) from the software can be used for informing clinical management for these criteria:

- Weight: No limitations in software
- Height: No limitations in software
- Health: No limitations in software
- Patient state: Patient is not user

## 2.4 Contraindications

Simpleware Medical is not intended for grading/sizing of anatomical defects to determine acceptable treatment plans. The output of the software is not a surgical plan itself, but helps with informing clinical management decisions – the software is not intended for use as a surgical plan.

No acceptance and performance testing, recurrent testing and maintenance is required for use. Software updates and patches will be provided by the manufacturer, where needed.

The DICOM Standard (version 3.0; 2022) is followed. Compliance details are available on the Simpleware website (<https://www.synopsys.com/simpleware/resources/regulatory.html>). We strongly recommend running through tests of the PACS export functionality on a PACS test system before using it in a clinical setting.

Simpleware Medical has the ability to process, store or discard information contained in medical image files such as DICOM files during the import process of these files. The import process can involve different data transfer methods including USB, CD/DVD, disk drive or server-based storage systems such as PACS. When such files contain personal patient information, it is the responsibility of the end-users to follow the local laws related to appropriate handling of personal data – for example HIPAA (USA) and GDPR (EU) – and to discard any information when required. Please refer

Symbol	Meaning
	<ul style="list-style-type: none"><li>• Hint or Information.</li><li>• Risk level: low.</li><li>• The provided information may be useful but ignoring it has a very low risk of causing patient injury or death.</li></ul>
	<ul style="list-style-type: none"><li>• Warning or potential usage error.</li><li>• Risk level: medium.</li><li>• The provided information is related to an operation that has succeeded but the user needs to be aware of details which may indicate undesirable side effects.</li></ul>

### 3.5.3 Recommended System Requirements

Operating System	Hardware
<ul style="list-style-type: none"> <li>• Windows 11</li> <li>• Windows Server 2019</li> </ul> <p>Only 64-bit operating systems are supported.</p>	<ul style="list-style-type: none"> <li>• Processor: Intel Core i7 or equivalent</li> <li>• Memory (RAM): 16 GB or higher**</li> <li>• Graphics card: Workstation grade OpenGL 3.2 compatible graphics card with 1024 MB VRAM</li> <li>• Screen resolution: 1920 x 1080; True color (32bit) or more</li> <li>• Disk space: 100 GB</li> </ul> <p>** As a guideline, allow twice as much RAM as your datasets size for volume rendering visualization and about 4 to 8 times for segmentation, surface mesh generation etc. The memory required depends heavily on the size and complexity of your image data and resultant model to be analyzed. Please contact Simpleware Support for further information with any specific queries.</p>

### 3.5.4 Uninstalling Simpleware Medical

Please refer to the steps for uninstalling programs on your particular version of Windows. In most cases, uninstalling can

## 4 Manufacturer Contact Details

### **Simpleware Medical V-2024.06**

Manufactured in June 2024 by Synopsys Simpleware Product Group

Synopsys (Northern Europe) Ltd.

Bradninch Hall, Castle Street

Exeter, Devon, EX4 3PL

United Kingdom

Phone: +44(0)1392 642635

Email: [simpleware@synopsys.com](mailto:simpleware@synopsys.com)

SRN: GB-MF-000002376

Simpleware Medical is a CE-marked product

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Synopsys (Nordeuropa) Ltd.
Synopsys International Ltd.
Synopsys Switzerland LLC